

em-tec Type	em-tec Part #	BioPharm Size	Tube Size [in inch]			Tube Size [in millimeter]		
			ID	WT	OD	ID	WT	OD
BCT 1/8 x 1/16"	12386	119 16	1/16 1/8	3/32 1/16	1/4 1/4	1.59 3.18	2.38 1.59	6.35 6.35
BCT 11/64 x 3/64"	12345	-	11/64*	3/64*	17/64*	4.30	1.25	6.80
BCT 3/16 x 1/16"	12344	120 25	1/8 3/16	3/32 1/16	5/16 5/16	3.18 4.76	2.38 1.59	7.94 7.94
BCT 1/4 x 1/16"	12341	15 17	3/16 1/4	3/32 1/16	3/8 3/8	4.76 6.35	2.38 1.59	9.53 9.53
BCT 1/4 x 3/32"	12342	24 18	1/4 5/16	3/32 1/16	7/16 7/16	6.35 7.94	2.38 1.59	11.11 11.11
BCT 3/8 x 1/16"	12346	26 -	1/4 3/8	1/8 1/16	1/2 1/2	6.35 9.53	3.18 1.59	12.70 12.70
BCT 3/8 x 3/32"	12367	36, 122	3/8	3/32	9/16	9.53	2.38	14.29
BCT 3/8 x 1/8"	12607	185 73	5/16 3/8	5/32 1/8	5/8 5/8	7.94 9.53	3.97 3.18	15.88 15.88
BCT 1/2 x 3/32"	12343	-	1/2	3/32	11/16	12.70	2.38	17.46
BCT 1/2 x 1/8"	12388	70, 86, 193 82	3/8 1/2	3/16 1/8	3/4 3/4	9.53 12.70	4.76 3.18	19.05 19.05
BCT 1/2 x 3/16"	12513	88 184	1/2 5/8	3/16 1/8	7/8 7/8	12.70 15.88	4.76 3.18	22.23 22.23
BCT 3/4 x 1/8"	12530	87 189 -	1/2 5/8 3/4	1/4 3/16 1/8	1 1 1	12.70 15.88 19.05	6.35 4.76 3.18	25.40 25.40 25.40
BCT 3/4 x 3/16"	12512	90, 191	3/4	3/16	1 - 1/8	19.05	4.76	28.56
BCT 1 x 1/8"	12613	91 -	3/4 1	1/4 1/8	1 - 1/4 1 - 1/4	19.05 25.40	6.35 3.18	31.75 31.75
BCT 1 x 3/16"	12558	92	1	3/16	1 - 3/8	25.40	4.76	34.93
BCT 1-1/4 x 3/16"	12698	-	1 - 1/4	3/16	1 - 5/8	31.75	4.76	41.28

*Approximate inch fractions



Where Innovation Flows

em-tec GmbH | Lerchenberg 20 | 86923 Finning, Germany | P: +49 8806 9236 0 | F: +49 8806 9236 50 | E: em-tec-info@psgdover.com | em-tec.de

em-tec BIOPROCESSING Tubing Guide V3.0 | EMT-10100-F-01-A4 © 2020 PSG, a Dover company

em-tec reserves the right to modify the information and illustrations contained in this document without prior notice. This is a not-contractual document.

