

# BioProTT™ FlowSU System

READY. CONNECT. FLOW.



## WE'VE LISTENED. WE'VE INNOVATED. WE'VE SOLVED YOUR PAIN POINTS.

Introducing the BioProTT™ FlowSU System – an innovative flow measurement system for single-use biopharma applications.

The new system features a first-of-its-kind disposable sensor that eliminates the need for calibration and a reusable and easy-to-mount holder profile with integrated electronics.

### READY

- Easy spec-in
- Easy order handling

### CONNECT

- Seamless integration
- Smart connection

### FLOW

- High level of features
- Low level of maintenance


## SMART & EASY FLOW MEASUREMENT INTEGRATION

### SINGLE-USE PART

- 1 • Standardized single-use sensor
  - Available in different sizes
  - Well-known ultrasonic technology inside
  - 2-in-1: Flow measurement & air-in-line detection with one sensor
  - Made of USP class VI material
- 2 • Easy integration into manifolds with standard connection
  - Straight flow path for seamless integration and less shear stress

### MULTI-USE PART

- 3 • Panel mounted re-usable holder profile
  - Integrated electronics and Power over Ethernet (PoE)
  - Cableless connection with the sensors
  - Status LED
- 4 • Toolless attachment with special locking technology
  - Station fits with every sensor size



#### READY

Easy Spec-In  
Easy Order Handling

#### CONNECT

Seamless Integration  
Smart Connection

#### FLOW

High Level of Features  
Low Level of Maintenance

PBT-E-30000-F-01

© 2023 PSG®, a Dover company

Authorized PSG® Partner:

### ProcessHQ, Inc.

410 N. Varney Street  
Burbank, CA 91502  
818 861-7641 (voice)  
818 475-1589 (fax)  
info@processHQ.com



PSG Biotech | em-tec GmbH  
Lerchenberg 20  
86923 Finning, Germany  
Tel: +49 8806 9236 0  
Fax: +49 8806 9236 50  
em-tec-sales@psgdover.com  
[psgdover.com/readyconnectflow](https://psgdover.com/readyconnectflow)

PSG® reserves the right to modify the information and illustrations contained in this document without prior notice. This is a non-contractual document.